



# recreo™



## Treatment - Simple, Fast and Safe

**STEP 1** - Press the on/ off button.

**STEP 2** - Use your hand to part your hair to expose your scalp. Place the **Recreo** on the scalp. The the bristles would help to seperate your hair easily.

**STEP 3** - Start at the front of your head move the **Recreo** device across your scalp.

**STEP 4** - When the **Recreo** beeps (every 4 seconds), pull the device away from your scalp and move onto another parting line.

**STEP 5** - Move the **Recreo** sideways.

**STEP 6** - Repeat step 1 to step 5 until you finish the treatment with 8-minute duration.

**DO NOT** Use the **Recreo** device more than 8 minutes at a time, 3 times a week.